

부서	대상 장비(기기코드)	공지 성능	공지 주기
지구환경 연구부	유도결합 플라즈마 시스템 (ICP System, AI11)	• Sensitivity : > 50 V /ppm • Accuracy : < $\pm 0.2\%$	연 2회 (1월, 7월)
	안정동위원소 질량분석기 (SIRMS, AI21)	• Standard Deviation: O : $\pm \leq 0.20\%$, H : $\pm \leq 1\%$	연 2회 (1월, 7월)
	환경방사능측정 시스템 (ERMS, AI31)	• Efficiency Variation : < 5 % (@352keV) • Resolution at 352 keV : < 1.3 keV	연 2회 (1월, 7월)
	열이온화 질량분석기 (TIMS, AI51)	• Abundance sensitivity < 2ppm at m/z 237 with respect to ^{238}U • External precision $^{87}\text{Sr}/^{86}\text{Sr}$ < 5ppm 1RSD of NIST SRM987	연 2회 (1월, 7월)
	다검출기 불활성기체 질량분석시스템 (SVMS, AI71)	• ^{40}Ar blank : < 8 fAmp • Mass resolution: > 200	연 2회 (1월, 7월)
	OSL 자동측정장비 (OSL, AI72)	• Beta source dose rate: Gy/s	연 4회 (3월, 6월, 9월, 12월)
	글로우방전 질량분석기 (GD-MS, AI74)	• Detection Limit: < 10ppb • Resolution: 10,000	연 2회 (1월, 7월)
	초고분해능 이차이온질량분석기 (SHRIMP, GT-01)	• Mass Resolution of $^{238}\text{U}^{16}\text{O}$ in high U-bearing zircon: > 5,000 • ^{206}Pb Sensitivity (cps/ppm/nA) in standard zircon of SL13: > 12	연 12회 (매월)
		• External precision for U-Pb zircon dating (using two standard zircons): < 2%	연 2회 (1월, 7월)
	레이저삭박 다검출기 유도결합플라즈마 질량분석기 (LA-MC-ICPMS, GT-03)	• Abundance sensitivity at mass 237: < 1ppm • Accuracy : external precision of standard zircons (91500, Plesovice) < 2%	연 2회 (1월, 7월)

분석지원장비 상태 보고서

Code (기기코드)	Instrument (장비명)
AI11	ICP System (유도결합 플라즈마 시스템)

Report Date (일시)	2018.07.03
------------------	------------

Sensitivity	Accuracy
Under repair	Under repair

분석지원장비 상태 보고서


Code (기기코드)	Instrument (장비명)
AI21	Stable Isotope Ratio Mass Spectrometer (안정동위원소 질량분석기)

Report Date (일시)	2018.07.03
------------------	------------

Standard Material (W-32593)


$\delta^2\text{H}_{\text{VSMOW}} = -55.2 \text{ ‰ (SD: } \leq \pm 1.0 \text{ ‰)}$

N	Sample	Corrected
1	W-32593	-55.8
2	W-32593	-55.4
3	W-32593	-55.7
4	W-32593	-56.4
5	W-32593	-55.3
6	W-32593	-54.6
7	W-32593	-54.6
8	W-32593	-55.3
9	W-32593	-55.2
10	W-32593	-54.4
	SD	0.62



$\delta^{18}\text{O}_{\text{VSMOW}} = -8.68 \text{ ‰ (SD: } \leq \pm 0.20 \text{ ‰)}$

N	Sample	Corrected
1	W-32593	-8.68
2	W-32593	-8.69
3	W-32593	-8.65
4	W-32593	-8.68
5	W-32593	-8.68
6	W-32593	-8.67
7	W-32593	-8.71
8	W-32593	-8.67
9	W-32593	-8.68
10	W-32593	-8.69
	SD	0.02




분석지원장비 상태 보고서


Code (기기코드)	Instrument (장비명)
AI31	Environmental Radioactivity Measurement System (환경방사능측정시스템)

Report Date (일시)	2018.07.03
------------------	------------

측정효율변화 (@352 keV)



분해능 (Resolution)



분석지원장비 상태 보고서

Code (기기코드)	Instrument (장비명)
AI51	Thermal Ionization Mass Spectrometer (열이온화질량분석기)

Report Date (일시)	2018.07.02
------------------	------------

Peak Plat profile
(<240ppm at H6 with respect to ^{82}Sr)

Accuracy
(Calibrator : External precision $^{87}\text{Sr}/^{86}\text{Sr} < 5\text{ppm}$
1 RSD of NIST SRM987)

Peak Flat Analysis



Normalised Results

Cycle: 21 / 20 Method: Peak H6
Baseline Type: Low & High Side

Collector: High6 (Faraday), 3000 steps
Mass: 81.331

Serial No.: 1310
Sample ID:

Sample No: 14




분석지원장비 상태 보고서

Code (기기코드)	Instrument (장비명)
AI71	Multi-collector static vacuum mass spectrometer (다검출기 비활성기체 질량분석기)

Report Date (일시)	2018.07.02
------------------	------------

Mass Resolution (at ^{40}Ar)



Mass resolution : 210

Blank

[illegible][illegible]



분석지원장비 상태 보고서

Code (기기코드)	Instrument (장비명)
AI74	Glow discharge mass spectrometer (고분해능 글로우방전 질량분석기)

Report Date (일시)	2018.07.02
------------------	------------

Sensitivity : 2.5×10^9 cps
Detection limits : < 0.5 ppb

Resolution($m/\Delta m$)
 $173.99887/0.015 = 11600$



분석지원장비 상태 보고서

Code (기기코드)	Instrument (장비명)
GT01	Sensitive High Resolution Ion MicroProbe (고분해능 이차이온 질량분석기)

Report Date (일시)	2018.07.03	²⁰⁶ Pb Sensitivity (cps/ppm/nA) in standard zircon of SL13: > 12	24.6
------------------	------------	---	------

Mass Resolution of ²³⁸ U ¹⁶ O in high U-bearing zircon: > 5,000	External precision for U-Pb zircon dating (using two standard zircons) : < 2%
Under repair	Under repair

분석지원장비 상태 보고서

Code (기기코드)	Instrument (장비명)
GT03	Laser Ablation Multi-Collector Inductively Coupled Plasma Mass Spectrometer (레이저삭박 다검출기 유도결합 플라즈마질량분석기)
Report Date (일시)	2018.07.02

Abundance sensitivity
(Calibrator : Abundance sensitivity <1ppm
at m/z 237 with respect to ^{238}U)

Accuracy
(Calibrator : External precision of standard
zircons(91500, Plesovice) < 2%)

Nu Plasma II Isotope Analysis Results

Date and Time: 2 July 2017 (15:05)
Sample Name: SRM112a-1(250ppb) abundance sensitivity with filters
Analysis Method File: C:\Nu Plasma II\Analyses\Standard Analyses\Abundance Sens
IC0 (U nat).nrf
Data File: C:\Nu Plasma II\Results\Data_2245.csv


Tau Correction invoked

Axial Mass(es) used: 245.317 248.490


Block Number: 1

	Mean	StdErr	Num Measurements
--	------	--------	------------------

Abundance Sensitivity (ppm)	8.709718E-01	3.58E-03	19 out of 20
234 U Intensity (Volts)	6.415416E-04	5.03E-06	18 out of 20
237 U Intensity (Volts)	7.104617E-06	3.28E-08	19 out of 20
234 U Intensity (Cps)	1.779161E+04	6.43E+01	18 out of 20
237 U Intensity (Cps)	6.871796E+02	2.05E+00	19 out of 20



standard 91500 zircon
(mean = 1062.4 ± 1.2)



standard Plesovice zircon
(concordia age = 337.48 ± 0.70)